

## Lab-Tips® - LTC125

### Large Closed Cell Foam Swab

100 ppi closed cell foam sock head swab

The LTC125 large rectangular swab is an excellent general-purpose swab for cleaning surfaces. Its closed cell polyurethane foam head possesses excellent durability and cushioning. A rigid internal paddle provides firm support.



#### Key Attributes

- Polyurethane foam over rigid paddle
- Closed cell structure
- Thermally bonded head
- Laundered and packaged in Class 10 {ISO Class 4} cleanroom

#### Benefits

- Low particles and fibres
- Good chemical resistance
- Excellent strength
- No contaminating adhesives
- Free from Silicone, amides and DOP

#### Applications

- Cleaning of sensitive equipment
- Applying lubricants and other liquids
- Removing excess materials
- Scrubbing applications

#### Physical Characteristics

<b>Head Material</b>	Closed Cell Foam	<b>Handle Thickness</b>	2.8mm {0.110"}
<b>Head Width</b>	13mm {0.511"}	<b>Handle Length</b>	100mm {3.94"}
<b>Head Thickness</b>	8mm {0.315"}	<b>Total Swab Length</b>	125mm {4.92"}
<b>Head Length</b>	25mm {0.984"}	<b>Head Bond</b>	Thermal
<b>Handle Material</b>	Polypropylene	<b>Handle Text</b>	Berkshire
<b>Handle Width</b>	5.5mm {0.216"}		

[www.berkshire.uk.com](http://www.berkshire.uk.com)

Contact: Tel + 44 1953 562800 / [info@berkshire.uk.com](mailto:info@berkshire.uk.com)

America Tel 1 413 528 2602

[info@berkshire.com](mailto:info@berkshire.com)

Europe Tel + 44 1953 562800

[enquiries@berkshire.uk.com](mailto:enquiries@berkshire.uk.com)

SE Asia Tel 65 6252 4313

[enquiries@berkshire.com.sg](mailto:enquiries@berkshire.com.sg)

Japan Tel 81 3 4530 9883

[master@berkshire.co.jp](mailto:master@berkshire.co.jp)

## Technical Data:

ATTRIBUTE		LTC125	TEST METHOD
Non-Volatile Residue	IPA mg/swab	1.4	IENT-RP-CC004.4 Sec 8.1.2
Ions (µg/swab)	Chloride	0.028	IENT-RP-CC004.4 Sec 8.2.2
	Sulfate	0.023	
	Nitrate	ND*	
	Phosphate	ND*	
	Fluoride	ND*	
	Potassium	0.014	
	Calcium	0.11	
	Sodium	0.0060	
	Magnesium	ND*	

## Notes:

- Data shown are typical values
- Due to differences in test methods applied and equipment utilised by different swab manufacturers, valid product comparisons may only be obtained through side-by-side testing in the same test facility, under similar conditions.
- \* ND\* – below the detection limit of the equipment.

## Order Information:

Product	Part Number	Size	Swabs/Pk	Pks/Cs	Case size
Large Closed Cell Foam	LTC125.5	125mm	100	5	220 x 160 X 105mm


## Other Berkshire Products



Wipers




Mop Systems



Face Masks



Glove Liners



Documentation Systems



Swabs