

# CONSUMABLES

UK BROCHURE :: 2023





## CONTENTS

INTRODUCTION	01	<b>CLOTHING</b>	<b>31</b>
REGIONAL ORDERING AND LABELS	02	DISPOSABLE COVERALL WITH INTEGRAL HOOD	32
<b>ALCOHOLS</b>	<b>03</b>	DISPOSABLE COVERALL WITH INTEGRAL HOOD AND FEET STERILE	33
HAND SANITISER	04	DISPOSABLE LABCOAT NON-STERILE	34
TRIGGER SPRAY	05	SINGLE USE CLEANROOM COMPATIBLE SOCK NON-STERILE	35
RFU 5L	06	DISPOSABLE FACEMASKS, BEARD SNOOD AND MOB CAP NON-STERILE	36
IPA POLYPROPYLENE TUB WIPE	07	<b>GLOVES</b>	<b>37</b>
IPA AND IMS POLYCELLULOSE POUCH WIPE	08	GLOVE DECISION TREE	38
IPA POLYESTER POUCH WIPE	09	POLYCHLOROPRENE 300mm HAND SPECIFIC GLOVE STERILE	39
IMS SILVER PACK POLYETHYLENE/POLYPROPYLENE PRE-SATURATED WIPE	10	LATEX 300mm HAND SPECIFIC GLOVE STERILE	40
<b>DISINFECTANTS</b>	<b>11</b>	NITRILE 300mm HAND SPECIFIC GLOVE STERILE	41
ALPHA, BETA AND DELTA TRIGGER SPRAY	12	NITRILE 405mm HAND SPECIFIC GLOVE STERILE	42
ALPHA, BETA, DELTA AND DI 5L RFU	13	NITRILE 230MM AMBIDEXTROUS GLOVE BLUE BOXED NON-STERILE	43
ALPHA, BETA AND DELTA POLYCELLULOSE POUCH WIPE	14	LATEX 300MM AMBIDEXTROUS GLOVE BULK PACKED NON-STERILE	44
<b>DISPOSABLE MOP SYSTEM</b>	<b>15</b>	NITRILE 300MM AMBIDEXTROUS GLOVE BULK PACKED NON-STERILE	45
ISOLATOR BONNET AND TOOL	16	NITRILE 300mm AMBIDEXTROUS GLOVE BULK PACKED NON-STERILE	46
IPA AND IMS DISPOSABLE MICROFIBRE MOP STERILE	17	<b>STATIONERY</b>	<b>47</b>
ALPHA, BETA AND DELTA DISPOSABLE MICROFIBRE MOP STERILE	18	CLEANROOM MARKER PEN BLACK AND BALLPOINT PEN STERILE	48
DRY DISPOSABLE MICROFIBRE MOP STERILE	19	CLEANROOM PAPER NON-STERILE	49
WASTE BAGS, MOP HEADS, TROLLEY AND OTHER ACCESSORIES	20	CLEANROOM MUNISING PAPER STERILE AND NON-STERILE	50
<b>DRY WIPES</b>	<b>21</b>	<b>OVERSHOES</b>	<b>51</b>
POLYCELLULOSE WIPE STERILE	22	POLYTHENE OVERSHOE BLUE NON-STERILE	52
POLYCELLULOSE WIPE NON-STERILE	23	PVC OVERSHOE NON-STERILE	53
POLYESTER WIPE STERILE	24	OVERSHOE ANTI-SLIP NON-STERILE	54
POLYESTER WIPE NON-STERILE	25	<b>MATS</b>	<b>55</b>
MICRODENIER WIPE NON-STERILE	26	ACCESS MAT MEDIUM TACK NON-STERILE	56
POLYCELLULOSE ROLL WIPE NON-STERILE	27	ACCESS MAT LOW TACK NON-STERILE	57
LAUNDERED DUAL LAYER POLYESTER WIPE STERILE	28	BLUE WALK ACCESS MAT WASHABLE NON-STERILE	58
<b>GOGGLES</b>	<b>29</b>	GREY PRO ACCESS MAT WASHABLE NON-STERILE	59
INDIVIDUALLY PACKED GOGGLE NON-STERILE	30	BLUE SELF-INSTALL ACCESS MAT NON-STERILE	60
		<b>RECYCLING INFORMATION</b>	<b>61</b>



---

## INTRODUCTION

---

“At Micronclean, our passion is to be the first to develop new technological solutions that change the shape of the markets we serve creating quality and efficiency for our customers.”

For over 40 years Micronclean has focused on contamination control products for companies in the pharma, healthcare, microelectronics, aerospace and semiconductor industries. As a result, we are incredibly proud of our cleanroom consumables range which we have sold to over 30 countries around the world in the last 5 years.

The 2023 edition of our Consumables Catalogue has been updated to include new products such as the Autoclavable Goggles, alongside recycling information for all our products.

A large percentage of our range is manufactured in the UK under cleanroom conditions. For those products we do not manufacture we have full control of the supply chain and set the specification. Having the capability to manufacture means we are able to offer bespoke products and maintain sufficient stock levels to meet customer demand and our customers are able to audit our facilities.

For more information on any of our product lines please visit our website [www.micronclean.com](http://www.micronclean.com) or contact us directly on [cs@micronclean.co.uk](mailto:cs@micronclean.co.uk) or +44 (0) 1754 767377

**Michelle Simpson**  
*Product Director*

## REGIONAL ORDERING AND LABELS

### MICRONCLEAN UPDATES ITS PRODUCT LABELS TO DEAL WITH DEMAND FROM EXPORT MARKETS.

Following increasing demand for our market leading disinfectant and cleaning products, we have had to enhance our existing product labels to ensure that we remain compliant.

As a result, we have generated an enhanced label which enables us to remain compliant with labeling regulations and allows us to carry up to three languages in addition to English.

So, the result of these changes is that the product code for our existing CleanGuard 1 IPA 500ml Trigger Spray Sterile will migrate from the current ZSAL64499 to ZSAL64499-A for region A orders, ZSAL64499-B for region B orders and so on. Note this regionalisation only affects products with active chemicals.

The new look labels with the associated regional codes for products will be rolled out progressively across our whole consumable range.

### Regional Labels in Use

For a UK consumable customer, the default option will always be that they order region A products and they will normally receive region A products. However, as the English language is present on all the labels and labelling regulations are common across region A-E, we might on occasion substitute the same product with a region B label for instance.

For a Micronclean customer in Italy for instance, the product should always be ordered for region B and to comply with the local language regulations we will supply the product with region B labelling.


Although English is present on all the labels, for regions G and I the products are destined for Australia and India respectively. This is because they have to conform to a different set of labeling regulations in each of these countries and as such these products can only be sold and used in these countries.

Please note if you have any questions please do not hesitate to contact us.

REGIONS	<b>A</b>	UK   Ireland   Germany   Austria
	<b>B</b>	UK   Spain   Portugal   Italy
	<b>C</b>	UK   France   Belgium   Holland
	<b>D</b>	UK   Sweden   Denmark   Norway
	<b>E</b>	UK   Greece   Poland
	<b>F</b>	UK   China
	<b>G</b>	UK   Australia
	<b>H</b>	Brazil
	<b>I</b>	UK   India
	<b>J</b>	UK   Croatia   Turkey

# ALCOHOLS





**RECYCLING INFORMATION**

+ info see page 62



## HAND SANITISER CLEANGUARD 1

**USAGE**

Designed as a sanitising hand rub for use in a professional environment

**BENEFITS**

- 70% alcohol content
- Manufactured with WFI
- UK BPR Registered
- Validated Performance:
  - EN 13727:2012+A2:2015 – Pass
  - EN 13624:2013 – Pass

**FEATURES**

- Transparent bottle
- Simple push down dispensing pump
- Fragrance-free
- Non-sticky
- Quick drying
- Supplied double bagged

Product Code	Description	Size	Pack Qty	Case Qty
ZNIP50025	IPA 500ml Hand Sanitiser Non-Sterile	500ml	Each	24 x 500ml



**RECYCLING INFORMATION**

+ info see page 62

## ALCOHOL TRIGGER SPRAY CLEANGUARD 1 & 2

### USAGE

Designed for use in cleanrooms  
Cleanroom transfer process  
Disinfection of hard surfaces  
Use within Isolators / LAFs

### BENEFITS

No bag in the bottle  
70% alcohol blended with WFI quality water  
Filtered to 0.2µm at point of fill  
Mechanically sealed trigger head  
Low centre of gravity bottle design  
Transparent bottle  
Long trigger handle, ease of use

### FEATURES

Trigger spray with spray / off feature  
2 year shelf life  
Gamma irradiation stable  
Supplied double bagged

Product Code	Description	Size	Pack Qty	Case Qty
ZSAL64499	IPA 500ml Trigger Spray Sterile	500ml	Each	24 x 500ml
ZSAL64950	IPA 950ml Trigger Spray Sterile	950ml	Each	15 x 950ml
ZNAL64950	IPA 950ml Trigger Spray Non-Sterile	950ml	Each	15 x 950ml
ZSAL65499	IMS 500ml Trigger Spray Sterile	500ml	Each	24 x 500ml
ZSAL65950	IMS 950ml Trigger Spray Sterile	950ml	Each	15 x 950ml
ZNAL65950	IMS 950ml Trigger Spray Non-Sterile	950ml	Each	15 x 950ml

The above product table does not show the regional product codes for clarity. However, to ensure that product labels with the correct languages are ordered and dispatched please ensure that the correct regional suffix is ordered as per the table on page 2.

**RECYCLING INFORMATION**

+ info see page 62



## IPA AND IMS 5 Litre RFU CLEANGUARD 1 & 2

**USAGE**

Disinfection of hard surfaces

**BENEFITS**

- 70% Alcohol blended with WFI quality water
- Bulk container for larger surface areas
- Easy to use with bucket mopping systems
- Filtered to 0.2µm at point of fill

**FEATURES**


- Tamper evident PP vented cap
- Gamma irradiation stable
- 2 year shelf life
- Supplied double bagged

Product Code	Description	Size	Pack Qty	Case Qty
ZSAL64005	IPA 5 Litre RFU Sterile	5L	Each	2 x 5L
ZSAL65005	IMS 5 Litre RFU Sterile	5L	Each	2 x 5L

The above product table does not show the regional product codes for clarity. However, to ensure that product labels with the correct languages are ordered and dispatched please ensure that the correct regional suffix is ordered as per the table on page 2.



RECYCLING INFORMATION



+ info see page 62



## IPA POLYPROPYLENE TUB WIPE VERIGUARD 1

### USAGE

Hard surface wiping particularly suited to outer areas

### BENEFITS

Designed for low grade cleanroom areas  
Designed for high volume usage  
Long lasting alcohol retention  
Re-fillable pouch available

### FEATURES

Continuous filament thermobonded nonwoven polypropylene wipe  
18gsm wipe weight  
2 years shelf life  
Supplied double bagged

Product Code	Description	Size	Pack Qty	Case Qty
WSVG05020	IPA Polypropylene Tub Wipe Sterile	200mm x 200mm	200	12 x 200
WSVG01020	IPA Polypropylene Tub Wipe Refill Sterile	200mm x 200mm	200	20 x 200
WNVG05020	IPA Polypropylene Tub Wipe Non-Sterile	200mm x 200mm	200	12 x 200

The above product table does not show the regional product codes for clarity. However, to ensure that product labels with the correct languages are ordered and dispatched please ensure that the correct regional suffix is ordered as per the table on page 2.



**RECYCLING INFORMATION**

+ info see page 62



## IPA AND IMS POLYCELLULOSE POUCH WIPE VERIGUARD 2

### USAGE

Hard surface wiping particularly suited to transfer processes

### BENEFITS

Increased alcohol retention from re-sealable pouch  
 High absorbency  
 Strong multifunctional wipe  
 Small pack size available for small batch manufacturing, large pack size suitable for large batch manufacturing

### FEATURES

White hydro-entangled Polycellulose wipe  
 68gsm wipe weight  
 2 years shelf life  
 Supplied double bagged

Product Code	Description	Size	Pack Qty	Case Qty
WSVG04023	IPA Polycellulose C-folded Pouch Wipes Sterile	230mm x 230mm	100	8 x 100
WSCW92009	IPA Sessional Polycellulose C-folded Pouch Wipes Sterile	230mm x 230mm	15	40 x 15
WNVG04023	IPA Polycellulose C-folded Pouch Wipes Non-Sterile	230mm x 230mm	100	8 x 100
WSVG09023	IMS Polycellulose C-folded Pouch Wipes Sterile	230mm x 230mm	100	8 x 100
WSCW91009	IMS Sessional Polycellulose C-folded Pouch Wipes Sterile	230mm x 230mm	15	40 x 15

The above product table does not show the regional product codes for clarity. However, to ensure that product labels with the correct languages are ordered and dispatched please ensure that the correct regional suffix is ordered as per the table on page 2.

**RECYCLING INFORMATION**

+ info see page 62



## IPA POLYESTER POUCH WIPE VERIGUARD 3

### USAGE

Cleanroom hard surface wipe  
Use in Isolators / Laminar flow cabinets

### BENEFITS

Increased alcohol retention from re-sealable pouch  
Laser cut edges  
Strong multifunctional wipe  
Small pack size available for small batch manufacturing, medium pack size suitable for larger batch manufacturing

### FEATURES

100% continuous filament polyester wipe  
2 years shelf life  
Supplied double bagged

Product Code	Description	Size	Pack Qty	Case Qty
WSVG06030	IPA Sessional Polyester V-folded Pouch Wipe Sterile	300mm x 300mm	15	30 x 15
MCZSLA50S	IPA Polyester Laundered Non-folded Pouch Wipes Sterile	300mm x 300mm	50	5 x 50
WNVG06030	IPA Sessional Polyester V-folded Pouch Wipe Non-Sterile	300mm x 300mm	15	30 x 15
MCZNLA50N	IPA Polyester Laundered Non-folded Pouch Wipes Non-Sterile	300mm x 300mm	50	5 x 50

The above product table does not show the regional product codes for clarity. However, to ensure that product labels with the correct languages are ordered and dispatched please ensure that the correct regional suffix is ordered as per the table on page 2.

**RECYCLING INFORMATION**

+ info see page 62



## IMS SILVER PACK POLYETHYLENE / POLYPROPYLENE PRE-SATURATED WIPE VERIGUARD 9

**USAGE**

Cleanroom hard surface wipe

**BENEFITS**

- Increased alcohol retention from re-sealable pouch
- Designed for large batch manufacturing
- Blended with de-ionised water (DI)
- Long lasting alcohol retention

**FEATURES**

- Nonwoven interleaved wipes
- 2 years shelf life
- Supplied double bagged

Product Code	Description	Size	Pack Qty	Case Qty
WSCW97ZZZ		200mm x 220mm	200	8 x 200

The above product table does not show the regional product codes for clarity. However, to ensure that product labels with the correct languages are ordered and dispatched please ensure that the correct regional suffix is ordered as per the table on page 2.

# DISINFECTANTS



**RECYCLING INFORMATION**

+ info see page 62



## ALPHA, BETA AND DELTA TRIGGER SPRAY CLEANGUARD 3, 4 AND 5

### USAGE

Disinfection and cleaning of hard surfaces

Annex 1 states that for disinfection to be effective it must be a two-stage process. The first stage should be a process to remove surface contamination and the second stage should be disinfection using a rotational system based on two types of disinfectants which have different modes of action

Micronclean has developed Delta as a sterile neutral detergent to remove surface contamination. This should be followed by Alpha and Beta as sterile rotational disinfectants which have different modes of action and prevent the long-term build-up of resistance

### BENEFITS

- No bag in the bottle
- Blended with WFI quality water
- Filtered to 0.2µm at point of fill
- Mechanically sealed trigger head
- Low centre of gravity bottle design
- Transparent bottle
- Long trigger handle, ease of use

### FEATURES

- Trigger spray with spray / off feature
- 2 year shelf life
- Gamma irradiation stable
- Supplied double bagged

Product Code	Description	Size	Pack Qty	Case Qty
ZSAP01500	Alpha 500ml Trigger Spray Sterile	500ml	Each	6 x 500ml
ZSBT02500	Beta 500ml Trigger Spray Sterile	500ml	Each	6 x 500ml
ZSND05950	Delta 950ml Trigger Spray Sterile	950ml	Each	15 x 950ml

The above product table does not show the regional product codes for clarity. However, to ensure that product labels with the correct languages are ordered and dispatched please ensure that the correct regional suffix is ordered as per the table on page 2.



RECYCLING INFORMATION



+ info see page 62

## ALPHA, BETA, DELTA AND DI 5L RFU CLEANGUARD 3, 4, 5 AND 6

### USAGE

Disinfection and cleaning of hard surfaces and floors

Annex 1 states that for disinfection to be effective it must be a two-stage process. The first stage should be a process to remove surface contamination and the second stage should be disinfection using a rotational system based on two types of disinfectants which have different modes of action

Micronclean has developed Delta as a sterile neutral detergent to remove surface contamination. This should be followed by Alpha and Beta as sterile rotational disinfectants which have different modes of action and prevent the long-term build-up of resistance

### BENEFITS

Blended with WFI quality water  
Bulk container for larger surface areas  
Easy to use with bucket mopping systems  
Filtered to 0.2µm at point of fill

### FEATURES

Tamper evident PP vented cap  
Gamma irradiation stable  
2 year shelf life  
Supplied double bagged

Product Code	Description	Size	Pack Qty	Case Qty
ZSAPO3005	Alpha 5 Litre RFU Sterile	5L	Each	2 x 5L
ZSBTO4005	Beta 5 Litre RFU Sterile	5L	Each	2 x 5L
ZSND05005	Delta 5 Litre RFU Sterile	5L	Each	2 x 5L
ZSDI82005	DI Water 5 Litre RFU Sterile	5L	Each	2 x 5L
ZNDI82005	DI Water 5 Litre RFU Non-sterile	5L	Each	2 x 5L

The above product table does not show the regional product codes for clarity. However, to ensure that product labels with the correct languages are ordered and dispatched please ensure that the correct regional suffix is ordered as per the table on page 2.



## ALPHA, BETA AND DELTA POLYCELLULOSE POUCH WIPE VERIGUARD 2

### USAGE

Hard surface wiping

Annex 1 states that for disinfection to be effective it must be a two-stage process. The first stage should be a process to remove surface contamination and the second stage should be disinfection using a rotational system based on two types of disinfectants which have different modes of action

Micronclean has developed Delta as a sterile neutral detergent to remove surface contamination. This should be followed by Alpha and Beta as sterile rotational disinfectants which have different modes of action and prevent the long-term build-up of resistance

### BENEFITS

- Increased liquid retention from re-sealable pouch
- Designed for small batch manufacturing
- High absorbency

### FEATURES


- White hydro-entangled Polycellulose wipe
- 68gsm wipe weight
- 2 years shelf life
- Supplied double bagged

Product Code	Description	Size	Pack Qty	Case Qty
WSAP01009	Alpha Sessional Polycellulose C-folded Pouch Wipe Sterile	230mm x 230mm	15	10 x 15
WSBT02009	Beta Sessional Polycellulose C-folded Pouch Wipe Sterile	230mm x 230mm	15	10 x 15
WSND05009	Delta Sessional Polycellulose C-folded Pouch Wipe Sterile	230mm x 230mm	15	10 x 15

The above product table does not show the regional product codes for clarity. However, to ensure that product labels with the correct languages are ordered and dispatched please ensure that the correct regional suffix is ordered as per the table on page 2.



# DISPOSABLE MOP SYSTEM





**RECYCLING INFORMATION**

+ info see page 62



## ISOLATOR BONNET AND TOOLS

### PUREGUARD 6 AND 7

**USAGE**

Designed specifically to help clean the inside of isolators, RABs and LAF cabinets

**BENEFITS**

Universal fit  
 Dual layer polyester guarantees high levels of absorbency  
 Lightweight isolator tool available in 2 handle lengths with 360° rotation

**FEATURES**

Low linting, dual-layer 100% continuous filament polyester for increased absorbency  
 250gsm bonnet weight  
 Supplied sterile one per bag ready for use  
 2 year sterile shelf life  
 Handle supplied Non-Sterile ready for sterilisation

Product Code	Description	Size	Pack Qty	Case Qty
ISED0140	Dry Dual Layer Polyester Laundered Isolator Bonnet Sterile	One Size	EA	40
ZNSS01035	Stainless Steel Isolator Cleaning Tool 35cm Non-Sterile	One size	EA	1
ZNSS01061	Stainless Steel Isolator Cleaning Tool 61cm Non-Sterile	One size	EA	1



**RECYCLING INFORMATION**

+ info see page 62

## IPA AND IMS DISPOSABLE MICROFIBRE MOP STERILE PUREGUARD 1 AND 2

### USAGE

Designed for use on cleanroom floors

### BENEFITS

Uniformly dosed to achieve a coverage of 15m<sup>2</sup>  
 No need for in cleanroom mixing or dosing  
 2 pocket design enables a hands-free and secure connection to the mop frame  
 Eliminates the risk associated with exposed chemicals and soiled water

### FEATURES

Proven 2 year shelf life  
 Individually packed mops are suitable for smaller cleanrooms  
 Multi packed mops are suitable for larger cleanrooms  
 Supplied double bagged

Product Code	Description	Size	Pack Qty	Case Qty
MSMI24016	IPA Individually Packed Disposable Microfibre Mop Sterile	440mm x 140mm	1	30 Packs of 1
MSMI21016	IPA Multi Packed Disposable Microfibre Mop Sterile	440mm x 140mm	10	3 Packs of 10
MSME24016	IMS Individually Packed Disposable Microfibre Mop Sterile	440mm x 140mm	1	30 Packs of 1
MSME21016	IMS Multi Packed Disposable Microfibre Mop Sterile	440mm x 140mm	10	3 Packs of 10

The above product table does not show the regional product codes for clarity. However, to ensure that product labels with the correct languages are ordered and dispatched please ensure that the correct regional suffix is ordered as per the table on page 2.



**RECYCLING INFORMATION**

+ info see page 62



## ALPHA, BETA AND DELTA DISPOSABLE MICROFIBRE MOP STERILE PUREGUARD 3, 4 AND 5

### USAGE

Annex 1 states that for disinfection to be effective it must be a two-stage process. The first stage should be a process to remove surface contamination and the second stage should be disinfection using a rotational system based on two types of disinfectants which have different modes of action

Micronclean has developed Delta as a sterile neutral detergent to remove surface contamination. This should be followed by Alpha and Beta as sterile rotational disinfectants which have different modes of action and prevent the long-term build-up of resistance

### BENEFITS

Uniformly dosed to achieve a coverage of 15m<sup>2</sup>

No need for in cleanroom mixing or dosing

2 pocket design enables a hands-free and secure connection to the mop frame

Eliminates the risk associated with exposed chemicals and soiled water

### FEATURES

Proven 2 year shelf life

Individually packed mops are suitable for smaller cleanrooms

Multi packed mops are suitable for larger cleanrooms

Supplied double bagged

Product Code	Description	Size	Pack Qty	Case Qty
MSMA24016	Alpha Individually Packed Disposable Microfibre Mop Sterile	440mm x 140mm	1	30 Packs of 1
MSMA21016	Alpha Multi Packed Disposable Microfibre Mop Sterile	440mm x 140mm	10	3 Packs of 10
MSMB24016	Beta Individually Packed Disposable Microfibre Mop Sterile	440mm x 140mm	1	30 packs of 1
MSMB21016	Beta Multi Packed Disposable Microfibre Mop Sterile	440mm x 140mm	10	3 Packs of 10
MSMN24016	Delta Individually Packed Disposable Microfibre Mop Sterile	440mm x 140mm	1	30 Packs of 1
MSMN21016	Delta Multi Packed Disposable Microfibre Mop Sterile	440mm x 140mm	10	3 Packs of 10

The above product table does not show the regional product codes for clarity. However, to ensure that product labels with the correct languages are ordered and dispatched please ensure that the correct regional suffix is ordered as per the table on page 2.



**RECYCLING INFORMATION**

+ info see page 62

## DRY DISPOSABLE MICROFIBRE MOP STERILE PUREGUARD 6

### USAGE

Designed for use on cleanroom floors

### BENEFITS

2 pocket design enables a hands-free and secure connection to the mop frame

### FEATURES

Proven 2 year shelf life  
 Individually packed mops are suitable for smaller cleanrooms  
 Multi packed mops are suitable for larger cleanrooms  
 Supplied double bagged

Product Code	Description	Size	Pack Qty	Case Qty
MSMD24016	Dry Individually Packed Disposable Microfibre Mop Sterile	One size	1	30 Packs of 1
MSMD21016	Dry Multi Packed Disposable Microfibre Mop Sterile	One Size	10	4 Packs of 10

The above product table does not show the regional product codes for clarity. However, to ensure that product labels with the correct languages are ordered and dispatched please ensure that the correct regional suffix is ordered as per the table on page 2.

Please see individual product specifications



RECYCLING INFORMATION



## WASTE BAGS, MOP HEADS, TROLLEY AND OTHER ACCESSORIES PUREGUARD 7 AND 8


### BENEFITS

Designed to work alongside the PureGuard Mopping System

Product Code	Description	Size	Pack Qty	Case Qty
BSP1C1012	Polyethylene Bag Sterile	One Size	10	12 Packs of 10
MTSS01001	Stainless Steel Standard Mop Trolley Non-Sterile	One size	Each	1
MFSS01001	Stainless Steel Lightweight Mop Frame 40cm Non-Sterile	One Size	Each	1
MFSS79226	Stainless Steel Shrouded Locking Mop Frame 40cm Non-Sterile	One Size	Each	1
MHSS02001	Stainless Steel Fixed Mop Handle 142cm Non-Sterile	One Size	Each	1
MHSS02002	Aluminium Adjustable Mop Handle 100cm-180cm Non-Sterile	One Size	Each	1
MTSS01001	Stainless Steel Standard Mop Trolley Non-Sterile	One Size	Each	1

The above product table does not show the regional product codes for clarity. However, to ensure that product labels with the correct languages are ordered and dispatched please ensure that the correct regional suffix is ordered as per the tables on page.

# DRY WIPES



**RECYCLING INFORMATION**

+ info see page 62



## POLYCELLULOSE WIPE STERILE WHITEGUARD 2

**USAGE**

Hard surface wipe  
Transfer of product into cleanrooms

**BENEFITS**

Excellent absorbency

**FEATURES**

Low linting, multipurpose wipes  
68gsm wipe weight  
Hydroentangled 55% Polyester, 45% Cellulose blend

Product Code	Description	Size	Pack Qty	Case Qty
WSWG04009	Polycellulose Non-Folded Wipe Sterile	230mm x 230mm (9" x 9")	100	10 Packs of 100
WSWG06009	Polycellulose Non-Folded Wipe Sterile	230mm x 230mm (9" x 9")	300	8 Packs of 300
WSWG04012	Polycellulose Non-Folded Wipe Sterile	300mm x 300mm (12" x 12")	100	10 Packs of 100
WSCS93009	Polycellulose Non-Folded Wipe Sterile	230mm x 230mm (9" x 9")	300	8 Packs of 300
WSCS93012	Polycellulose Non-Folded Wipe Sterile	300mm x 300mm (12" x 12")	10 x 10	17 Packs of 100



**RECYCLING INFORMATION**

+ info see page 62



## POLYCELLULOSE WIPE NON-STERILE WHITEGUARD 2

### USAGE

Hard surface wipe  
Transfer of product into cleanrooms


### BENEFITS

Excellent absorbency

### FEATURES

Low linting, multipurpose wipes  
68gsm wipe weight  
Hydroentangled 55% Polyester, 45% Cellulose blend

Product Code	Description	Size	Pack Qty	Case Qty
WNWG04009	Polycellulose Non-Folded Wipe Non-Sterile	230mm x 230mm (9" x 9")	100	10 Packs of 100
WNWG06009	Polycellulose Non-Folded Wipe Non-Sterile	230mm x 230mm (9" x 9")	300	16 Packs of 300
WNWG04012	Polycellulose Non-Folded Wipe Non-Sterile	300mm x 300mm (12" x 12")	100	20 Packs of 100
WNCS93009	Polycellulose Non-Folded Wipe Non-Sterile	230mm x 230mm (9" x 9")	300	8 Packs of 300
WNCS93012	Polycellulose Non-Folded Wipe Non-Sterile	300mm x 300mm (12" x 12")	100	17 Packs of 100



+ info see page 63



## POLYESTER WIPE STERILE WHITEGUARD 3

### USAGE

Hard surface wipe  
Transfer of product into cleanrooms

### BENEFITS

Ultra low linting wipe

### FEATURES

Laser sealed edges to prevent linting  
100% knit continuous filament polyester  
Medium levels of absorbency

Product Code	Description	Size	Pack Qty	Case Qty
WSTG05009	Polyester Non-Folded Wipe Sterile	230mm x 230mm (9" x 9")	150	10 packs of 150
WSTG03012	Polyester Non-Folded Wipe Sterile	300mm x 300mm (12" x 12")	50	12 packs of 50
WSPW80009	Polyester Laundered Non-Folded Wipe Sterile	230mm x 230mm (9"x9")	150	10 packs of 150
MCDELN50S	Polyester Laundered Single Layer Non-Folded Wipe Sterile	300mm x 300mm (12"x12")	50	10 packs of 50

RECYCLING INFORMATION



+ info see page 63



## POLYESTER WIPE NON-STERILE WHITEGUARD 3

### USAGE

Hard surface wipe  
Transfer of product into cleanrooms

### BENEFITS

Ultra low linting wipe

### FEATURES

Laser sealed edges to prevent linting  
100% knit continuous filament polyester  
Medium levels of absorbency

Product Code	Description	Size	Pack Qty	Case Qty
WNTG04009	Polyester Non-Folded Wipe Non-Sterile	230mm x 230mm (9" x 9")	100	10 packs of 100
WNTG05009	Polyester Non-Folded Wipe Non-Sterile	230mm x 230mm (9" x 9")	150	10 packs of 150
WNTG05012	Polyester Non-Folded Wipe Non-Sterile	300mm x 300mm (12" x 12")	150	10 packs of 150
WNTG03012	Polyester Non-Folded Wipe Non-Sterile	300mm x 300mm (12" x 12")	50	10 packs of 50


**RECYCLING INFORMATION**  
 + info see page 63



## MICRODENIER WIPE NON-STERILE WHITEGUARD 5

**USAGE**

Hard surface wipe  
Excellent wipe for lens and screen cleaning

**BENEFITS**

Ultra low linting wipe

**FEATURES**

100% continuous filament nylon/polyester wipe

Product Code	Description	Size	Pack Qty	Case Qty
WNMD04009	Microdenier Non-folded Wipe Non-Sterile	230mm x 230mm (9" x 9")	100	10 Packs of 100



RECYCLING INFORMATION



+ info see page 62

## POLYCELLULOSE ROLL WIPE NON-STERILE WHITEGUARD 7

### USAGE

For use in support areas  
Hard surface wipe


### BENEFITS

Excellent absorbency

### FEATURES

Low linting, multipurpose wipes  
Log roll of 400 wipes

Product Code	Description	Size	Pack Qty	Case Qty
WNPW81010	Polycellulose Roll Wipe Non-Sterile White	254mm x 381mm	400 sheets/roll	2 x 1 Roll
WNPB81010	Polycellulose Roll Wipe Non-Sterile Blue	254mm x 381mm	400 sheets/roll	2 x 1 Roll


**RECYCLING INFORMATION**

+ info see page 63



## LAUNDERED DUAL LAYER POLYESTER WIPE STERILE WHITEGUARD 11

### USAGE

Hard surface wipe  
Transfer of product into cleanrooms


### BENEFITS

### FEATURES

Dual layer 100% continuous filament polyester  
Laser sealed edges to prevent linting  
Medium levels of absorbency

Product Code	Description	Size	Pack Qty	Case Qty
MHPELN50S	Dual Layer Polyester Laundered Non-folded Wipe Sterile	300mm x 300mm (12" x 12")	50	4 Packs of 50

# GOGGLES



RECYCLING INFORMATION

+ info see page 63



## INDIVIDUALLY PACKED GOGGLE NON-STERILE VISIONGUARD 1

### USAGE

Specifically designed for use in Cleanroom Environments

### BENEFITS

Fully compatible with the requirements of Annex 1

### FEATURES

Supplied in a Non-Sterile condition but packed in an Autoclave bag ready for immediate sterilisation

40 process lifespan for goggles without significant fogging

Lightweight, supplied with adjustable head strap and easy-to-use snap-on clip

Indirect ventilation to top of goggle to prevent fogging and protect from splashes

Compatible with most prescription glasses


EN 166 3 4 BT Compliant

CE and UKCA Marked

Product Code	Description	Size	Pack Qty	Case Qty
ONACL0112	Individually Packed Goggle Non-Sterile	One Size	1	



# CLOTHING





RECYCLING INFORMATION

+ info see page 63



## DISPOSABLE COVERALL WITH INTEGRAL HOOD

### SUREGUARD 1

#### USAGE

Designed for cleanroom use

#### BENEFITS

CE Marked to PPE Category III Type 5 & 6  
Anti-static

#### FEATURES

Single use cleanroom coverall  
CE and UKCA Marked.  
Polypropylene with polyethylene laminate fabric  
Integrated hood  
Latex free elastic at hood, wrists and ankles  
Available in Sterile and Non-Sterile Forms

Product Code	Description	Size	Pack Qty	Case Qty
DSWH04	Disposable Coverall with Integral Hood Sterile	SML - 2XL	1	25
DNWH04	Disposable Coverall with Integral Hood Non-Sterile	SML - 2XL	1	25



RECYCLING INFORMATION



+ info see page 63

## DISPOSABLE COVERALL WITH INTEGRAL HOOD AND FEET STERILE SUREGUARD 3

### USAGE

Designed for cleanroom use

### BENEFITS

CE Marked to PPE Category III Type 5 & 6  
Anti-static

### FEATURES

Single use cleanroom coverall  
CE and UKCA Marked.  
Polypropylene with polyethylene laminate fabric  
Integrated hood  
Latex free elastic at hood, wrists and ankles

Product Code	Description	Size	Pack Qty	Case Qty
DSWH21	Disposable Coverall with Integral Hood and Feet Sterile	SML - 2XL	1	25



**RECYCLING INFORMATION**

+ info see page 63




---

## DISPOSABLE LABCOAT NON-STERILE SUREGUARD 2

---

**USAGE**

Designed for laboratory and cleanroom use

**BENEFITS**

CE Marked to PPE Category III Type 5 & 6  
Anti-static







**FEATURES**

Single use cleanroom labcoat  
CE and UKCA Marked  
Polypropylene with polyethylene laminate fabric  
Latex free elastic at wrist

Product Code	Description	Size	Pack Qty	Case Qty
DNWH10	Disposable Labcoat Non-Sterile	SML - 2XL	1	25



**RECYCLING INFORMATION**

-  
-  
-  

+ info see page 63

## SINGLE USE CLEANROOM COMPATIBLE SOCK NON-STERILE TUFFGUARD 1

**FEATURES**

Specifically developed for use in cleanrooms to help with compliance with Annex 1

**BENEFITS**

Sock shaped for wearer comfort

**FEATURES**

White sock shaped  
100% synthetic fibres  
One Size

Product Code	Description	Size	Pack Qty	Case Qty
SNWH50001	Single Use Cleanroom Compatible Sock Non-Sterile	One Size	100 (50 pairs)	100 (50 pairs)

**RECYCLING INFORMATION**

+ info see page 63



**FACEGUARD 1**  
FACEMASK WITH TIES

**FACEGUARD 2**  
FACEMASK WITH EARLOOPS

**FACEGUARD 3**  
FACEMASK WITH HEADLOOPS

**FACEGUARD 5 AND 6**  
BEARD SNOOD AND MOP CAP

**DISPOSABLE FACEMASK FACEGUARD 1, 2 AND 3**  
**BEARD SNOOD AND MOB CAP NON-STERILE FACEGUARD 5 AND 6**

**FEATURES**

DISPOSABLE FACEMASKS

- One Size
- Available with ties, earloops or headloops
- Available in Sterile and Non-Sterile forms

BEARD SNOOD AND MOP CAP


- One Size
- Comfort fit

Product Code	Description	Size	Pack Qty	Case Qty
MNTW01100	Multi Packed Tie-On Face Mask Non-Sterile	One Size	100	10 x 100
MNTW01020	Individually Packed Tie-On Face Mask Non-Sterile	One Size	20	50 x 20
MSTW01020	Individually Packed Tie-On Face Mask Sterile	One Size	20	15 x 20
MNEW01100	Multi Packed Earloop Face Mask Non-Sterile	One Size	100	10 x 100
MNEW01020	Individually Packed Earloop Face Mask Non-Sterile	One Size	20	50 x 20
MSHW01020	Individually Packed Earloop Face Mask Sterile	One Size	50	20 x 50
MSEW01020	Individually Packed Headloop Face Mask Sterile	One Size	20	25 x 20
MNHW01100	Multi Packed Headloop Face Mask Non-Sterile	One Size	100	10 x 100
MNHW01020	Individually Packed Headloop Face Mask Non-Sterile	One Size	20	50 x 20
MNWH01ZZZ	Beard Snood White Non-Sterile	One Size	100	5 X 100
HNWH01ZZZ	Mob Cap White Non-Sterile	One Size	100	10 X 100


# GLOVES



# GLOVE DECISION TREE



Do you want a **Sterile** or **Non-Sterile** Glove?







**RECYCLING INFORMATION**

+ info see page 63

## POLYCHLOROPRENE 300mm HAND SPECIFIC GLOVE PACKED IN L/R WALLET STERILE SKINGUARD 1

### USAGE

Specifically developed for use in cleanrooms to provide enhanced protection from cytotoxic drugs without compromising dexterity

Compatible with ISO Class 4+

### BENEFITS

Cytotoxic protection  
Enhanced Dexterity  
Latex free

### FEATURES

100% Polychloroprene  
CE / UKCA Marked  
Packed in L/R Wallet  
Hand specific  
Powder free  
Textured Palm  
Beaded cuff  
5 year shelf life

Product Code	Description	Size	Pack Qty	Case Qty
GSPG31	Polychloroprene 300mm Hand Specific Glove Packed in L/R Wallet Sterile	6.0 - 9.0	10 Pairs	5 Packs of 10 Pairs



RECYCLING INFORMATION

+ info see page 63



## LATEX 300mm HAND SPECIFIC GLOVE PACKED IN L/R WALLET STERILE SKINGUARD 5

### USAGE

Specifically developed for use in cleanrooms to provide enhanced dexterity, flexibility and wearer comfort

Compatible with ISO Class 4+

### BENEFITS

Enhanced Dexterity

### FEATURES

- 100% Latex
- CE / UKCA Marked
- Packed in L/R Wallet
- Hand specific
- Powder free
- Textured Palm
- Beaded cuff
- 5 year shelf life

Product Code	Description	Size	Pack Qty	Case Qty
GSLN05	Latex 300mm Hand Specific Glove Packed in L/R Wallet Sterile	6.0 - 9.0	10 Pairs	



**RECYCLING INFORMATION**

+ info see page 63

## NITRILE 300mm HAND SPECIFIC GLOVE PACKED IN L/R WALLET STERILE SKINGUARD 6

**USAGE**

Specifically developed for use in cleanrooms to provide dexterity, flexibility and wearer comfort  
Compatible with ISO Class 4+

**BENEFITS**

Good Dexterity  
Latex Free

**FEATURES**

CE / UKCA Marked  
Packed in L/R Wallet  
Hand specific  
Powder free  
Textured Palm  
Beaded cuff  
5 year shelf life

Product Code	Description	Size	Pack Qty	Case Qty
GSNW06	Nitrile 300mm Hand Specific Glove Packed in L/R Wallet Sterile	6.0 - 9.0	10 Pairs	20 Packs of 10 Pairs



**RECYCLING INFORMATION**

+ info see page 63



## NITRILE 405mm HAND SPECIFIC GLOVE PACKED IN L/R WALLET STERILE SKINGUARD 7

**USAGE**

Specifically developed for use in cleanrooms to provide dexterity, flexibility and wearer comfort  
Compatible with ISO Class 4+

**BENEFITS**

Good Dexterity  
Latex Free

**FEATURES**

CE / UKCA Marked  
Packed in L/R Wallet  
Hand specific  
Powder free  
Textured Palm  
Beaded cuff  
5 year shelf life

Product Code	Description	Size	Pack Qty	Case Qty
GSNW07	Nitrile 405mm Hand Specific Glove Packed in L/R Wallet Sterile	6.0 - 9.0	10 Pairs	20 Packs of 10 Pairs



**RECYCLING INFORMATION**

+ info see page 63

## NITRILE 230mm AMBIDEXTROUS GLOVE BLUE BOXED NON-STERILE SKINGUARD 2

**USAGE**

Specifically developed for use in cleanrooms  
Compatible with ISO Class 7/8


**BENEFITS**

Latex Free  
Bulk packed in cardboard dispensers

**FEATURES**

Blue  
Powder free  
Nitrile material  
Ambidextrous  
Textured surface for increased grip  
Beaded cuff

Product Code	Description	Size	Pack Qty	Case Qty
GNNB02	Nitrile 230mm Ambidextrous Glove Blue Boxed Non-Sterile	SML - XLG	150 Gloves	10 Packs of 150



**RECYCLING INFORMATION**

+ info see page 63



## LATEX 300mm AMBIDEXTROUS GLOVE BULK PACKED NON-STERILE SKINGUARD 10

**USAGE**

Specifically developed for use in cleanrooms to provide enhanced dexterity, flexibility and wearer comfort  
Compatible with ISO Class 4+

**BENEFITS**

Enhanced Dexterity

**FEATURES**

- 100% Latex
- CE / UKCA Marked
- Powder free
- Textured Palm
- Beaded cuff
- Bulk Packed

Product Code	Description	Size	Pack Qty	Case Qty
GNLN10	Latex 300mm Ambidextrous Glove Bulk Packed Non-Sterile	SML - XLG	100 Gloves	10 Packs of 100



**RECYCLING INFORMATION**

+ info see page 63

## NITRILE 300mm AMBIDEXTROUS GLOVE BULK PACKED NON-STERILE SKINGUARD 11

**USAGE**

Specifically developed for use in cleanrooms to provide dexterity, flexibility and wearer comfort  
Compatible with ISO Class 4+


**BENEFITS**

Good Dexterity  
Latex Free  
Liquid particle count <600 at 0.5µm

**FEATURES**

CE / UKCA Marked  
Packed in L/R Wallet  
Hand specific  
Powder free  
Textured Palm  
Beaded cuff  
Bulk Packed

Product Code	Description	Size	Pack Qty	Case Qty
GNNW11	Nitrile 300mm Ambidextrous Glove Bulk Packed Non-Sterile	SML - XLG	100 Gloves	10 Packs of 100



**RECYCLING INFORMATION**

+ info see page 63



## NITRILE 300mm AMBIDEXTROUS GLOVE BULK PACKED NON-STERILE SKINGUARD 12

**USAGE**

Specifically developed for use in cleanrooms to provide dexterity, flexibility and wearer comfort  
Compatible with ISO Class 6+

**BENEFITS**

Good Dexterity  
Latex Free  
Liquid particle count <3000 at 0.5µm


**FEATURES**


CE / UKCA Marked  
Packed in L/R Wallet  
Hand specific  
Powder free  
Textured Palm  
Beaded cuff  
Bulk Packed

Product Code	Description	Size	Pack Qty	Case Qty
GNNW12	Nitrile 300mm Ambidextrous Glove Bulk Packed Non-Sterile	SML- XLG	100 Gloves	10 Packs of 100



# STATIONERY





**RECYCLING INFORMATION**

+ info see page 63



**GUARDMARK 1**  
CLEANROOM BALLPOINT PEN STERILE



**GUARDMARK 2**  
CLEANROOM MARKER PEN BLACK STERILE

## CLEANROOM BALLPOINT PEN STERILE GUARDMARK 1 CLEANROOM MARKER PEN BLACK STERILE GUARDMARK 2

**USAGE**

**CLEANROOM BALLPOINT PEN**

Designed for controlled environments, this sterile ballpoint pen features low smear, low sodium black ink and inherent ESD properties

**CLEANROOM MARKER PEN**

Proven ability to write on most typical laboratory surfaces including plastic, glass, and stainless steel

**BENEFITS**

**CLEANROOM BALLPOINT PEN**

Low smear black ink

**CLEANROOM MARKER PEN**

Low particulate

**FEATURES**

**CLEANROOM BALLPOINT PEN**

Low smear black ink with low sodium content  
Inherent ESD properties  
Packed in ISO 14644-1 ISO Class 6 cleanroom  
2 years shelf life

**CLEANROOM MARKER PEN**

Quick-drying ink which is water and UV-light resistant  
Does not dry out even if the cap is left off for a day  
Packed in ISO 14644-1 ISO Class 6 cleanroom  
2 years shelf life

Product Code	Description	Size	Pack Qty	Case Qty
ZSBK11210	Cleanroom Ballpoint Pen Sterile	One Size	12	120
ZSBK03010	Cleanroom Marker Pen Black Sterile	One Size	1	10



**RECYCLING INFORMATION**

+ info see page 64

## CLEANROOM PAPER NON-STERILE GUARDMARK 3

### USAGE

Designed for in cleanroom use with laser or inkjet printers


### BENEFITS

Available in 5 colours: White, Blue, Pink, Green and Yellow  
Low particulate generation

### FEATURES

250 sheets per pack  
Basic weight 72gsm

Product Code	Description	Size	Pack Qty	Case Qty
ZNWH04072	Cleanroom Paper White Non-Sterile	A4	250 sheets	10 x 250
ZNBL04072	Cleanroom Paper Blue Non-Sterile	A4	250 sheets	10 x 250
ZNPI04072	Cleanroom Paper Pink Non-Sterile	A4	250 sheets	10 x 250
ZNGR04072	Cleanroom Paper Green Non-Sterile	A4	250 sheets	10 x 250
ZNYE04072	Cleanroom Paper Yellow Non-Sterile	A4	250 sheets	10 x 250



## CLEANROOM MUNISING PAPER STERILE AND NON-STERILE GUARDMARK 4

### USAGE

Designed for in cleanroom use with laser or inkjet printers

### BENEFITS


Low particulate generation  
Good heat stability  
ROHS Compliant


### FEATURES

Sterile and Non-Sterile versions available in a variety of colour (see table below)  
Manufactured from a combination of Cellulose Fibre and Synthetic Latex  
250 sheets per pack  
Basic weight 80gsm

Product Code	Description	Size	Pack Qty	Case Qty
ZNWH04022	Cleanroom Munising Paper White Non-Sterile	A4	250 sheets	10 x 250
ZNBL04022	Cleanroom Munising Paper Blue Non-Sterile	A4	250 sheets	10 x 250
ZNYE04022	Cleanroom Munising Paper Yellow Non-Sterile	A4	250 sheets	10 x 250
ZSWH04225	Cleanroom Munising Paper White Sterile	A4	250 sheets	12 x 250
ZSGR04022	Cleanroom Munising Paper Green Sterile	A4	250 sheets	12 x 250
ZSBL04225	Cleanroom Munising Paper Blue Sterile	A4	250 sheets	12 x 250

# OVERSHOES





+ info see page 64



## POLYTHENE OVERSHOE BLUE NON-STERILE TUFFGUARD 2

### FEATURES


Elastic opening for easy donning

Available in two sizes

Available in two materials depending on usage

Manufactured from PE

Product Code	Description	Size	Pack Qty	Case Qty
SNBL01016	Blue Heavy Duty Polythene Overshoe 16" Non-Sterile	420mm x 160mm	100	10 Packs of 100
SNBL31014	Blue Polythene Overshoe 14" Non-Sterile	360mm x 150mm	100	20 Packs of 100
SNBL31016	Blue Polythene Overshoe 16" Non-Sterile	410mm x 150mm	100	20 Packs of 100


**RECYCLING INFORMATION**



+ info see page 64

## PVC OVERSHOE NON-STERILE TUFFGUARD 3

### FEATURES

- Non slip sole
- Available in two sizes
- Manufactured from PVC
- Elastic opening for easy donning

Product Code	Description	Size	Pack Qty	Case Qty
SNWH81006	PVC 6"x14" Overshoe White Non-Sterile	6" x 14"	500	500
SNWH83016	PVC 6"x16" Overshoe White Non-Sterile	6" x 16"	500	500
SNWH82009	PVC 9"x14" Overshoe White Non-Sterile	9" x 14"	500	500
SNWH81016	PVC 9"x16" Overshoe White Non-Sterile	9" x 16"	500	500
SNRE81009	PVC 9"x16" Overshoe Red Non-Sterile	9" x 16"	500	500



+ info see page 64



## OVERSHOE ANTI-SLIP NON-STERILE

### TUFFGUARD 6

#### BENEFITS

Anti slip sole


#### FEATURES


Elastic opening for easy donning  
 Manufactured from Polypropylene (+EVA)  
 laminate

Product Code	Description	Size	Pack Qty	Case Qty
SNWH03016	Overshoe White Anti-Slip Non-Sterile	180mm x 400mm	100	10 Packs of 100



# MATS





RECYCLING INFORMATION

+ info see page 64



## ACCESS MAT MEDIUM TACK NON-STERILE GUARDTACK 1

### USAGE

Peel-Off access mat designed to capture contamination before entry into a controlled area

### BENEFITS

Medium "Tack" formulation gives a higher level of dirt capture than GuardTack 2 whilst not inhibiting free movement

### FEATURES

30 Individually numbered LDPE peel-off sheets per contamination control mat  
8 mats per pack  
Available in three colours (Blue / Grey / White)  
Available in three sizes

Product Code	Description	Size	Pack Qty	Case Qty
ANXX08018	Contamination Control Mats Medium Tack Non-Sterile	458mm x 1145mm (18" x 45")	30 sheets/mat	8 Mats
ANXX08026	Contamination Control Mats Medium Tack Non-Sterile	660mm x 1145mm (26" x 45")	30 sheets/mat	8 Mats
ANXX08336	Contamination Control Mats Medium Tack Non-Sterile	914mm x 914mm (36" x 36")	30 sheets/mat	8 Mats

Where XX = BL for Blue, GY for Grey and WH for White

RECYCLING INFORMATION



+ info see page 64



## ACCESS MAT LOW TACK NON-STERILE GUARDTACK 2

### USAGE

Peel-Off access mat designed to capture contamination before entry into a controlled area

### BENEFITS


Low "Tack" formulation gives a lower level of dirt capture compared to GuardTack 1 but does not affect free movement

### FEATURES

30 Individually numbered LDPE peel-off sheets per contamination control mat  
4 mats per pack  
Available in three colours (Blue / Grey / White)  
Available in five sizes


Product Code	Description	Size	Pack Qty	Case Qty
ANXX04018	Contamination Control Mats Low Tack Non-Sterile	458mm x 1145mm (18" x 45")	30 sheets/mat	4 Mats
ANXX04236	Contamination Control Mats Low Tack Non-Sterile	660mm x 914mm (24" x 36")	30 sheets/mat	4 Mats
ANXX04026	Contamination Control Mats Low Tack Non-Sterile	660mm x 1145mm (26" x 45")	30 sheets/mat	4 Mats
ANXX04036	Contamination Control Mats Low Tack Non-Sterile	914mm x 1145mm (36" x 45")	30 sheets/mat	4 Mats
ANXX04336	Contamination Control Mats Low Tack Non-Sterile	914mm x 914mm (36" x 36")	30 sheets/mat	4 Mats

Where XX = BL for Blue, GY for Grey and WH for White



**RECYCLING INFORMATION**

+ info see page 64



## BLUE WALK ACCESS MAT WASHABLE NON-STERILE GUARDTACK 3

**USAGE**

Washable, high performance contamination control mat for use with medium traffic loads

**BENEFITS**

Specially designed and manufactured to capture, retain and reduce foot and wheel-borne viable and non-viable particulates, over an extended period of time

**FEATURES**


Installed by trained service technicians.  
 3-5 year in service life  
 Silver based antimicrobial embedded into the surface  
 Max 1000kg per mat area

Product Code	Description	Size	Pack Qty	Case Qty
ANBLWLK12	Blue Walk Access Mat 1.2m Non-Sterile	1200mm	Each	1 Mat
ANBLWLK20	Blue Walk Access Mat 2.0m Non-Sterile	1200mm	Each	1 Mat

RECYCLING INFORMATION



+ info see page 64



## GREY PRO ACCESS MAT WASHABLE NON-STERILE GUARDTACK 4

### USAGE

Washable, high performance contamination control mat for use with high traffic loads


### BENEFITS

Specially designed and manufactured to capture, retain and reduce foot and wheel-borne viable and non-viable particulates, over an extended period of time

### FEATURES


Installed by trained service technicians  
 3-5 year in service life  
 Silver based antimicrobial embedded into the surface  
 Max 3000kg per mat area  
 Available with stainless steel edging to help protect the entry and exit off the mat surface

Product Code	Description	Size	Pack Qty	Case Qty
ANGYPRO12	Grey Pro Access Mat 1.2m Non-Sterile	1200mm	Each	1 Mat
ANGYPRO20	Grey Pro Access Mat 2.0m Non-Sterile	2000mm	Each	1 Mat



**RECYCLING INFORMATION**

+ info see page 64



## BLUE SELF-INSTALL ACCESS MAT NON-STERILE GUARDTACK 5

### USAGE

Washable, high performance contamination control mat for use with light traffic loads

### BENEFITS


Specially designed and manufactured to capture, retain and reduce foot and wheel-borne viable and non-viable particulates, over an extended period of time

### FEATURES

Available in 3 sizes  
 Self-install  
 3-5 year in service life  
 Silver based antimicrobial embedded into the surface  
 Max 400kg per mat area








Product Code	Description	Size	Pack Qty	Case Qty
ANBLINS06	Blue Self-install Access Mat 0.6m x 1.2m Non-Sterile	1200mm x 600mm	Each	1 Mat
ANBLINS10	Blue Self-install Access Mat 1.0m x 1.2m Non-Sterile	1000mm x 1200mm	Each	1 Mat
ANBLINS15	Blue Self-install Access Mat 1.0m x 1.5m Non-Sterile	1000mm x 1500mm	Each	1 Mat

# RECYCLING INFORMATION











# RECYCLING INFORMATION









**HAND SANITISER**  
Page 04

Component Name	Component Material	Disposal
 Bottle Pump	Mixed Materials	 Non-recyclable
 500ml Bottle	PET	 Recyclable
 Polybag	LDPE	 Recyclable
 Box	Cardboard	 Recyclable











**TRIGGER SPRAY**  
Page 05 | 12

Component Name	Component Material	Disposal
 Trigger Spray Head with Dip Tube	Mixed Materials	 Non-recyclable
 500ml/950ml Bottle	HDPE	 Recyclable
 Polybag	LDPE	 Recyclable
 Box	Cardboard	 Recyclable









**5 LITRE RFU**  
Page 06 | 13

Component Name	Component Material	Disposal
 Container Cap	HDPE	 Recyclable
 5L Container	PP	 Recyclable
 Polybag	LDPE	 Recyclable
 Box	Cardboard	 Recyclable









**POLYPROPYLENE TUB WIPES**  
Page 07

Component Name	Component Material	Disposal
 Wipe	Polypropylene	 Recyclable
 Tub	HDPE	 Recyclable
 Tub Lid	PE	 Recyclable
 Polybag	LDPE	 Recyclable
 Box	Cardboard	 Recyclable









**POLYCELLULOSE POUCH WIPES**  
Page 08 | 10 | 14

Component Name	Component Material	Disposal
 Wipe	Polycellulose	 Non-recyclable
 Wipe Pouch	Polyester and Polyethylene blend	 Non-recyclable
 Polybag	LDPE	 Recyclable
 Box	Cardboard	 Recyclable









**IPA POLYESTER POUCH WIPES**  
Page 09

Component Name	Component Material	Disposal
 Wipe	Polyester	 Recyclable
 Wipe Pouch	Polyester and Polyethylene blend	 Non-recyclable
 Polybag	LDPE	 Recyclable
 Box	Cardboard	 Recyclable









**IMS SILVER PACK POLYETHYLENE/POLYPROPYLENE WIPES**  
Page 10

Component Name	Component Material	Disposal
 Wipe	Polyethylene / Polypropylene	 Non-recyclable
 Wipe Pouch	OPP + PE + VMPET	 Non-recyclable
 Polybag	LDPE	 Recyclable
 Box	Cardboard	 Recyclable







**ISOLATOR BONNET AND TOOLS**  
Page 16

Component Name	Component Material	Disposal
 Bonnet	Polyester	 Recyclable
 Handle	Stainless Steel	 Recyclable
 Polybag	LDPE	 Recyclable
 Box	Cardboard	 Recyclable

**IPA AND IMS DISPOSABLE MICROFIBRE MOP STERILE**  
Page 17 | 18 | 19


Component Name	Component Material	Disposal
 Mop	Microfibre (Polyester and Polyamide)	 Non-recyclable
 Mop Pouch (Individually Packed Only)	Polyester and Polyethylene blend	 Non-recyclable
 Polybag	LDPE	 Recyclable
 Box	Cardboard	 Recyclable







**POLYCELLULOSE WIPES / ROLL WIPES**  
Page 22 | 23 | 27


Component Name	Component Material	Disposal
 Wipe	Hydroentangled 55% Polyester, 45% Cellulose blend	 Non-recyclable
 Polybag	LDPE	 Recyclable
 Box	Cardboard	 Recyclable












# RECYCLING INFORMATION







 **POLYESTER WIPES STERILE**  
Page 24 | 25 | 28

Component Name	Component Material	Disposal
 Wipe	100% Polyester	 Recyclable
 Polybag	LDPE	 Recyclable
 Box	Cardboard	 Recyclable

 **DISPOSABLE FACEMASKS, BEARD SNOOD AND MOB CAP**  
Page 36


Component Name	Component Material	Disposal
 Mobcap/Beard Snood	PP Fabric With Rubber Elastic	 Non-recyclable
 Facemask		 Non-recyclable
 Polybag	LDPE	 Recyclable
 Box	Cardboard	 Recyclable









 **MICRODENIER WIPES NON-STERILE**  
Page 26


Component Name	Component Material	Disposal
 Wipe	100% continuous filament nylon/ polyester wipe	 Non-recyclable
 Polybag	LDPE	 Recyclable
 Box	Cardboard	 Recyclable









 **POLYCHLOROPRENE GLOVES**  
Page 39


Component Name	Component Material	Disposal
 Gloves	Polychloroprene	 Recyclable
 L & R Wallet	HDPE	 Recyclable
 Polybag	HDPE	 Recyclable
 Box	Cardboard	 Recyclable







 **INDIVIDUALLY PACKED GOGGLE NON-STERILE**  
Page 30


Component Name	Component Material	Disposal
 Goggle	Thermoplastic Rubber, Silicone, Polycarbonate	 Non-recyclable
 Autoclave Bag	French Paper With Blue laminated Film (PE and CPP blend)	 Non-recyclable
 Outer Bag	Nylon	 Recyclable
 Box	Cardboard	 Recyclable









 **LATEX GLOVES**  
Page 40 | 44


Component Name	Component Material	Disposal
 Gloves	Latex	 Recyclable
 L & R Wallet	HDPE	 Recyclable
 Polybag	MDPE	 Recyclable
 Box	Cardboard	 Recyclable







 **DISPOSABLE COVERALL / DISPOSABLE LAB COAT**  
Page 32 | 33 | 34


Component Name	Component Material	Disposal
 Coverall	White Microporous PP with PE laminate with zip and rubber elastic	 Non-recyclable
 Polybag	LDPE	 Recyclable
 Box	Cardboard	 Recyclable







 **NITRILE GLOVES**  
Page 41 | 42 | 43 | 45 | 46

Component Name	Component Material	Disposal
 Gloves	Nitrile Butadiene Rubber	 Recyclable
 L & R Wallet	HDPE	 Recyclable
 Polybag	LDPE	 Recyclable
 Box	Cardboard	 Recyclable


 **CLEANROOM POLYESTER SOCKS**  
Page 35







Component Name	Component Material	Disposal
 Socks	Polyester/ HDPE	 Recyclable
 Polybag	LDPE	 Recyclable
 Box	Cardboard	 Recyclable


 **CLEANROOM MARKER PEN BLACK**  
Page 48







Component Name	Component Material	Disposal
 Pen	Polypropylene / Polyester	 Non-recyclable
 Polybag	LDPE	 Recyclable
 Box	Cardboard	 Recyclable


# RECYCLING INFORMATION







 **CLEANROOM PAPER**  
Page 49


Component Name	Component Material	Disposal
 Paper	Paper	 Recyclable
 Polybag	LDPE	 Recyclable
 Box	Cardboard	 Recyclable







 **PVC OVERSHOE / OVERSHOE ANTI-SLIP**  
Page 53 | 54


Component Name	Component Material	Disposal
 Overshoe	PVC With Rubber Elastic Band	 Non-recyclable
 Polybag	LDPE	 Recyclable
 Box	Cardboard	 Recyclable









 **CLEANROOM MUNISING PAPER**  
Page 50


Component Name	Component Material	Disposal
 Paper	Cellulose Fiber and Synthetic Latex	 Non-recyclable
 Polybag	LDPE	 Recyclable
 Box	Cardboard	 Recyclable







 **ACCESS MAT MEDIUM AND LOW TACK**  
Page 56 | 57

Component Name	Component Material	Disposal
 Mat	LDPE coated with water-based acrylic adhesive	 Non-recyclable
 Polybag	LDPE	 Recyclable
 Box	Cardboard	 Recyclable

 **POLYTHENE OVERSHOE BLUE**  
Page 52

Component Name	Component Material	Disposal
 Overshoe (Heavy Duty)	Chlorinated PE (CPE) With A Rubber Elastic Band	 Non-recyclable
 Overshoe	LDPE/ LLDPE with a Rubber Elastic Band	 Non-recyclable
 Polybag	LDPE	 Recyclable
 Box	Cardboard	 Recyclable

 **ACCESS MAT**  
Page 58 | 59 | 60

Component Name	Component Material	Disposal
 Mat		 Non-recyclable
 Polybag		 Recyclable
 Box		 Recyclable

**Recycling information relates to**

- UK Law and Practices
- Clean items that are free from contamination and residues

---

## NOTES

---



Roman Bank | Skegness | Lincolnshire | PE25 1SQ | UK  
T: +44 (0)1754 767 377  
E: sales@micronclean.co.uk  
W: www.micronclean.com



MCB 08/2024